



THANK YOU

For additional support & information, visit www.redux.com



RE DUX

#530-0021 v1.0
Instruction for use for 900-VD11 / DB11L-BL



RE DUX
Veri-Dry™

www.redux.com


USER MANUAL

Welcome to the Veri-Dry™ experience!

Veri-Dry™ is designed to protect what matters most: your devices. Our innovative, heat-free drying system uses advanced Electronic Rice™ technology to safely remove moisture from hearing aids, earbuds, and other delicate electronics.

We combine smart sensors, reusable materials, and precision engineering to make moisture removal simple, fast, and reliable.

Please read this manual carefully before first use and keep it for future reference. Your feedback helps us improve every detail. For support or comments, contact us at support@redux.com.



Protect. Reactivate. Reuse.

— The Redux™ Team —

DESIGNED AND ENGINEERED IN THE USA

1

Product Overview

The Veri-Dry™ System is made up of the Veri-Dry™ Case and DryBUDDY™. Veri-Dry™ is a portable, heat-free moisture-removal system designed to rescue and protect small electronic devices such as smartphones, hearing aids, earbuds, cochlear implants, and other hearables. It uses Electronic Rice™ desiccant beads to gently absorb moisture without heat. The built-in monitoring system tracks drying progress. It also indicates when reactivation is required, ensuring reliable, reusable performance. Keep DryBUDDY™ sealed in its Veri-Dry™ Case when not in use so it does not continue absorbing ambient moisture.

Warnings & Precautions

⚠ WARNING

- **INGESTION HAZARD:** This product contains a button cell or coin battery.
- **DEATH** or serious injury can occur if ingested.
- A swallowed button cell or coin battery can cause **Internal Chemical Burns** in as little as **2 hours**.
- **KEEP** new and used batteries **OUT OF REACH** of **CHILDREN**
- **Seek immediate medical attention** if a battery is suspected to be swallowed or inserted inside any part of the body.



*For more information go to page 16

Safety Information

Read all safety information before use. Failure to follow these instructions may result in injury or damage.

2

Electrical & Fire Safety

- Do not immerse the Veri-Dry™ case, DryBUDDY™, the USB-C cable, or reactivation sleeve in water.
- Do not disassemble; There are no user-serviceable parts.
- Only use the USB cable that came with your device and a 5V, minimum 10W output USB block.

Heat & Handling Safety









- The reactivation sleeve and Electronic Rice™ may become hot during reactivation. Allow to cool before handling. Wait 10 minutes after reactivation is complete before touching or moving the unit.
- Make sure the DryBUDDY™ moisture absorption slots are not blocked or completely covered when placed inside the Veri-Dry™ case.
- If the moisture absorption slots are covered, the device cannot effectively remove moisture.

Electromagnetic Considerations

- The reactivation sleeve contains a small magnet. While it poses no known electromagnetic interference risk, avoid placing near sensitive medical or magnetic equipment (e.g., pacemakers, MRI machines) as a precaution.

Keep for Reference

- Retain this manual for future use. Do not discard. It includes detailed instructions and troubleshooting information.



3

4

Product Components

Component	Label	Description
Power Switch	A	When looking at the front of DryBUDDY™ (and the green button is facing you), the off position is to the left and on is to the right.
Green Start Button	B	Press to begin the monitoring process when the Power Switch is ON.
LED Indicator	C	Dual-color LED (green/red) displays drying and reactivation status. See LED Indicator Guide, page 7.
Electronic Rice™ Chamber	D	Contains the moisture-absorbing Electronic Rice™ beads. These beads are sealed within the chamber and can be reactivated when saturated.
USB-C Port	E	External power input used exclusively for reactivating the Electronic Rice™ beads.
Battery Compartment	F	Holds one CR2032 coin cell battery secured behind a screw-lock cover. The battery powers the DryBUDDY™'s monitoring function.

5

Moisture Absorption Slots

G

Vents located on the bottom of the device. Do not block or cover these openings during operation.

Reactivation Sleeve

H

Special sleeve that fits over the DryBUDDY™ during reactivation to restore the Electronic Rice™ beads.

VeriDry™ Case

I

The case in which DryBUDDY™ stored and put your devices to dry.

Case Slide Latch

J

Slide latch to the left to close and to the right to open.

DryBUDDY™

K

DryBUDDY™ unit.

What's in the Box

- Veri-Dry™ Case
- DryBUDDY™ Unit
- Reactivation Sleeve
- USB-C Cable
- USB-C to USB-A Adapter
- User Manual
- Quick Start Guide

OFF ← → ON



6

LED Indicator Guide

State	Green LED	Red LED	Meaning / Action
Monitoring	Blinking	OFF	DryBUDDY™ is actively monitoring moisture levels inside the Veri-Dry™ Case. A typical drying cycle lasts up to 3-6 hours, depending on the amount of moisture present. The unit automatically goes to sleep off after 20 hours.
Dry Complete	Solid ON	OFF	The monitored devices are dry. The green LED remains ON for several hours after completion, then turns off automatically. Remove your devices and slide the Power Switch to OFF as the DryBUDDY™ to save battery life.
Sleep	OFF	OFF	If the switch is left ON and both lights are off, the DryBUDDY™ has gone to sleep. Slide the Power Switch to OFF to fully power down and conserve the CR2032 battery.

Reactivation Required

OFF

Solid ON

The Electronic Rice™ beads are saturated and cannot absorb any more moisture. Reactivation is required before the next use. Follow the Reactivation Procedure on page 9.

Reactivation (in sleeve)

OFF

Blinking

The Electronic Rice™ beads are reactivating inside the Reactivation Sleeve. This process takes approximately 3 hours.

Reactivation Complete

Solid ON

OFF

The reactivation process is complete. The Electronic Rice™ beads are fully refreshed and ready for use. Allow 10 minutes for cooling before removing the DryBUDDY™.

Sensor Fault

Solid ON

Solid ON

Turn power switch OFF, wait 10 seconds, then power ON again. If the issue persists, contact Redux™ Support at support@redux.com.

Initial Setup

Power Source

- **Battery:** A CR2032 coin cell battery comes pre-installed. It powers the DryBUDDY™ monitoring function during use inside the Veri-Dry™ Case.
- **USB-C Power:** Plug the supplied USB-C cable into the port (E) and a standard USB 5V adapter or USB port. USB-C power does not recharge the CR2032 battery, it is only used for reactivation.

Electronic Rice™

When new or freshly reactivated, the Electronic Rice™ beads appear bright orange. As they absorb moisture, they gradually darken indicating it's time for reactivation.

Starting a Dry Cycle

Follow these steps carefully to ensure optimal drying performance.

Power On DryBUDDY™

- Slide the Power Switch to ON, then press the Green Start Button to activate the monitoring process before placing it inside the Veri-Dry™ Case. DryBUDDY™ uses Electronic Rice™ and will always absorb moisture. It is important to keep DryBUDDY™ in its Veri-Dry™ Case, or other airtight packaging.

Insert Devices

- Place your electronic devices such as earbuds (e.g., AirPods), or hearing aids inside the Veri-Dry™ Case along with the powered-on DryBUDDY™ unit.
- Ensure your devices are powered OFF before closing Veri-Dry™ Case.



9

- Ensure devices do not interfere with the Veri-Dry™ Case seal. Close the Veri-Dry™ case lid. Slide latch to the left to close.

Drying Process

- Typical drying time is 3–6 hours, depending on the size of the device and the amount of moisture present. When drying is complete, the green LED will turn solid and remain lit for several hours before automatically turning off as the DryBUDDY™ goes to sleep. After drying, turn the Power Switch OFF and keep DryBUDDY™ inside the sealed Veri-Dry™ Case to prevent unnecessary moisture absorption.

Reactivation Needed

- If a solid red LED appears, it means the Electronic Rice™ beads are saturated and cannot absorb any more moisture.
- Follow the Reactivation Procedure (see page 11) to restore full drying capacity before next use.

Completing the Dry Cycle

Automatic Stop

- DryBUDDY™ stops monitoring automatically once moisture removal is complete—or after 20 hours, whichever comes first.
- A typical drying cycle lasts 3–6 hours, depending on humidity and device moisture levels.

Dry Complete Indicator

- When drying is complete, the green LED turns solid for several hours, then turns off automatically. This indicates that DryBUDDY™ no longer detects moisture and your devices are safe to use.
- After use, slide the Power Switch to OFF on the DryBUDDY™ to conserve battery life.



10

Remove Devices

- Once the green LED is solid or turns off, remove your devices from the Veri-Dry™ Case.

Storage Tip

- After every use, reseal DryBUDDY™ inside the Veri-Dry™ Case to prevent it from absorbing ambient moisture.

Electronic Rice™ Health

Reactivation Required Alert


- If the red LED turns solid ON during operation or at startup, the Electronic Rice™ beads are saturated and must be reactivated.
- The beads will also appear dark orange or brown, indicating full moisture absorption.
- To restore performance, turn the Power Switch OFF and follow the Reactivation Procedure section before your next use.

When to Reactivate

- Always reactivate before your next drying cycle to ensure full moisture-absorption performance.

How to Reactivate

- Slide the Power Switch to OFF.
- Plug the supplied USB-C cable into the port and connect it to a standard 5V USB adapter.
- Insert DryBUDDY™ fully into the Reactivation Sleeve.



- The Red LED will blink during the reactivating process.
- Keep DryBUDDY™ connected for approximately 3 hours while the reactivation process completes.
- When the green LED turns ON, reactivation is complete.
- Allow DryBUDDY™ to cool for at least 10 minutes before handling or using again.
- Store DryBUDDY™ in the sealed Veri-Dry™ Case when not in use to prevent unwanted moisture absorption.



Troubleshooting

Problem	Possible Cause	Solution
Device won't power on	Dead/missing CR2032	Replace battery and restart.
Red LED solid	Beads saturated	Perform Reactivation Procedure.
Both LEDs solid	Sensor fault	Power OFF, wait 10 seconds, restart.

Cleaning & Maintenance

- Turn OFF and disconnect USB-C before cleaning.
- Wipe exterior with a dry cloth.
- Do not use liquids or abrasive materials.
- Store DryBUDDY™ sealed in Veri-Dry™ Case.

Battery Handling & Replacement

The DryBUDDY™ is powered by a single CR2032 coin cell battery, which supports the monitoring function during use. The battery is not rechargeable.

When to Replace the Battery

- Replace the battery if DryBUDDY™ does not power ON or the LEDs fail to illuminate when switched ON.
- Always turn OFF the unit when not in use to conserve battery life.

How to Replace the Battery

- **Turn Power OFF:** Slide the Power Switch to OFF before replacing the battery.
- **Locate the Battery Compartment:** The battery compartment is located on the front of the DryBUDDY™.
- **Remove Screw:** Use a small Phillips screwdriver to remove the upper-right screw securing the battery door. Do not remove any other screws.
- **Replace Battery:** Remove the old CR2032 coin cell and insert a new one, positive (+) side facing up.
- **Secure Battery Door:** Reattach the door and tighten the screw gently. Do not overtighten.
- **Confirm the battery compartment is fully closed and secured.** If it does not close securely, stop using the product, remove the battery, and keep it away from children.
- **Test the Unit:** Slide the Power Switch to ON to confirm the LEDs light up correctly.

Disposal & Recycling

- Do not dispose of DryBUDDY™ or its accessories with household waste.
- Take the DryBUDDY™, Veri-Dry™ Cases reactivation sleeve, and USB-C cable to an approved e-waste recycling center for proper disposal.

Battery

- Remove the CR2032 coin cell battery before disposing of DryBUDDY™.
- Dispose of used batteries separately at a certified battery collection or recycling facility, following local regulations.

Warranty & Service

Your Veri-Dry™ system is covered by a 12-month limited warranty from the date of purchase.

Coverage

Covers defects in materials and workmanship under normal consumer use.

Exclusions

Does not cover damage resulting from misuse, unauthorized repairs, disassembly, or use of non-approved accessories.

How to Obtain Support

1. Keep proof of purchase (such as a receipt).
2. Contact Redux™ Support at support@redux.com for assistance.
3. Follow the instructions provided by the Support Team for warranty or service requests.

Condition

Specification

Operating Temperature

Operates between 5° C and 40° C (41° F to 104° F).

Reactivation Temperature

Internal components of the DryBUDDY™ may reach 220° F (104° C) during reactivation. The exterior surface of the reactivation sleeve may reach up to 140 degrees F (60 degrees Celsius) but remains safe to handle when inside the insulated sleeve.

Storage Temperature

Store between -20° C and 60° C (-4° F to 140° F) in a dry environment.

Battery Specifications

Battery Type: CR2032 (coin cell)
Nominal Battery Voltage: 3 VDC
Battery is non-rechargeable.

Operating Humidity

Functions effectively in non-condensing environments with 0-90% relative humidity.

WARNING: INGESTION HAZARD COIN/BUTTON CELL BATTERY

- Remove and immediately recycle or dispose of used batteries according to local regulations and keep away from children. Do NOT dispose of batteries in household trash or incinerate.
- Even used batteries may cause severe injury or death.
- Call a local poison control center for treatment information. For treatment information, call the National Battery Ingestion Hotline at 1-800-498-8666.
- Non-rechargeable batteries are not to be charged.
- Do not force discharge, recharge, disassemble, heat above 212°F (100°C), or incinerate. Doing so may result in injury due to venting, leakage or explosion resulting in chemical burns.
- Ensure the batteries are installed correctly according to polarity (+ and -).
- Do not mix old and new batteries, different brands or types of batteries, such as alkaline, carbon-zinc, or rechargeable batteries.
- Always completely secure the battery compartment. If the battery compartment does not close securely, stop using the product, remove the batteries, and keep them away from children.
- This product contains silica gel desiccant beads. Choking Hazard DO NOT EAT. Keep desiccant materials away from children.

Manufacturer Information

Revive Electronics, LLC dba Redux™
9715 Kincaid Drive, Suite 1200
Fishers, IN 46037 USA
Phone: 1-844-REDUX-IT
Email: support@redux.com
Web: www.redux.com/support

Patent Info:

*U.S. Patent Nos. 12,215,925; 12,276,454; 12,281,847.
Additional U.S. & international patents pending.*

